

# K2 RAPID TEST

## ONE-STEP DRUG SCREEN

The One-Step Synthetic Cannabinoid Drug Screen Test is a lateral flow immunoassay for the specific, qualitative detection of synthetic cannabinoid metabolites in human urine.



### ONE-STEP K2 RAPID TEST

Compound	Conc. (ng/ml)	Compound	Conc. (ng/ml)
JWH-018 Pentanoic acid	30	JWH-073 Butanoic acid	15
JWH-018-N-4-hydroxypentyl	200	JWH-073-N-4- hydroxybutyl	300
JWH-081-N-5- hydroxypentyl	1000	JWH-200-N-6-hydroxyindole	300
AM-2201-N-4-hydroxypentyl	1000	JWH-250-N-5-hydroxyindole	300
RCS-4-N-5- Carboxypentyl	250		

RESULTS YOU CAN TRUST

**POINT OF CARE MEDICAL**  
DIAGNOSTIC TESTING PRODUCTS

**877.836.5713** [NOBLEMEDICAL.COM](http://NOBLEMEDICAL.COM)  
19525 Janacek Court, Ste. 104 | Brookfield, WI 53045

