

6-MAM RAPID TESTS

RESULTS YOU CAN TRUST



WHY TEST FOR 6-MAM?

- 6-Monoacetylmorphine (6-MAM) is one of three active metabolites of Heroin (diacetylmorphine)
- 6-MAM is available in urine for no more than 24 hours and is created from Heroin.
- Thus, 6-MAM is the most specific marker for Heroin.

WHAT IS SPECIAL ABOUT IMPROVED 6-MAM TEST FROM NOBLE?

- Most of the tests for heroin are affected by Morphine, which raises concerns about use of 6-MAM Tests for testing Heroin abuse.
- Improved Noble 6-MAM Tests have undergone stringent validation testing and show no cross-reactivity with Morphine at 100 µg/ml, Oxycodone at 100 µg/ml, Codeine at 100 mg/ml and Hydrocodone at 100 µg/ml.
- Thus, Improved Noble 6-MAM Tests provide a highly accurate testing solution for Heroin abuse.

RESULTS YOU CAN TRUST

POINT OF CARE MEDICAL
DIAGNOSTIC TESTING PRODUCTS

877.836.5713 NOBLEMEDICAL.COM
19525 Janacek Court, Ste. 104 | Brookfield, WI 53045

