

ETG TEST

ETHYL GLUCURONIDE (ETG) RAPID URINE TEST

MORE ACCURATE. Strong indicator of alcohol ingestion within the previous 3-4 days.

PREFECT FOR ALCOHOL TREATMENT PROGRAMS.

Can give early warning of potential relapse.

Useful Window of
Detection of
8 to 80 Hours



Parameter: ETG



Calibrator: Ethyl Glucuronide

INTERPRETATION OF RESULTS

C
T

POSITIVE: Only one colored band appears, in the control region (C). No apparent colored band appears in the test region (T).

C
T

NEGATIVE: Two colored bands appear on the membrane. One band appears in the control region (C) and another band appears in the test region (T).

c
T

INVALID: Control band fails to appear. Results from any test which has not produced a control band at the specified read time must be discarded. Please review the procedure and repeat with a new test. If the problem persists, discontinue using the kit immediately and contact your local distributor.



RESULTS YOU CAN TRUST

POINT OF CARE MEDICAL
DIAGNOSTIC TESTING PRODUCTS

877.836.5713 NOBLEMEDICAL.COM
19525 Janacek Court, Ste. 104 | Brookfield, WI 53045

NOBLE
MEDICAL INC